

# *ACTIVATE*

PLACENTAL- DERIVED ALLOGRAFT

**NIVISION**™  
BIOMEDICAL TECHNOLOGIES

Exclusively Distributed by  
**WVC**  
Innovative Wound Consultants



# ACTIVATE

Activate™ Placental-Derived Allograft is comprised of the amniotic and chorionic membranes, and includes the intermediate (spongy) layer of the placenta. This product is designed as an innovative solution to serve as a barrier or wound cover across a spectrum of clinical needs. Leveraging the structural and biochemical properties of placental tissue, Activate™ offers a comprehensive biologic approach for a variety of acute or chronic wound applications.

## PLACENTAL TISSUE

Placental tissue, previously considered biologic waste following delivery, is a pristine newly formed extracellular matrix (ECM) to protect and support fetal development. Therefore, this tissue represents an untainted extracellular matrix.....essentially a neomatrix.

## DESCRIPTIONS OF LAYERS

Broadly, placental tissue includes an internal or fetus facing amnion and an external or maternal facing chorion. The physical and biochemical properties of the two primary membranes are subtly different owing to the position they occupy in the placenta. The amnion is a thinner yet more compact layer of extracellular matrix while the chorion is a slightly thicker but less compact matrix. Both layers are constituted of collagen fibers (predominantly Types I, III, IV and V), along with other ECM proteins (including laminin, fibronectin, proteoglycans and hyaluronic acid).

## SAFETY

Donated tissues for Activate™ products are collected from fully consented mothers undergoing full term c-sections. Each donor is screened according to the strict standards required by the U.S. Food and Drug Administration and the American Association of Tissue Banks. Further, Activate™ Placental-Derived Allograft products undergo a validated terminal sterilization process to help ensure these products are safe. All donated tissue is obtained in partnership with FDA regulated and accredited recovery organizations. Additional details of screening procedures can be found in the product package insert or Instructions for Use.





### PRODUCT SUMMARY

Activate™ Placental-Derived Allograft products are minimally manipulated and dehydrated and are derived from human amniotic and chorionic membranes and includes the intermediate (spongy) layer of the placenta. Activate™ allograft tissues retain the structural and functional characteristics of the starting placental membranes. The final products are packaged in different sizes and verified to be terminally sterilized to a  $10^{-6}$  SAL.



### INSTRUCTION FOR USE

Activate™ Placental-Derived Allograft (361 HCT/Ps) are intended as a natural biologic wound covering or skin substitute for cutaneous wounds. Use of Activate™ Placental-Derived Allograft by qualified health care professionals is for application in a physician office, outpatient, or inpatient setting.



### ORDERING INFORMATION

Activate™ tissue is available in the following sizes for a variety of choices depending on the patient and circumstance: Product HCPCS Level II code Q4301, Activate™ Matrix, per square centimeter.

#### PRODUCT SKU:

#### PRODUCT SIZE:

**ACT-27222**

**2 X 2**

**ACT-27224**

**2 X 4**

**ACT-27244**

**4 X 4**

**ACT-27248**

**4 X 8**

### REFERENCES:

- 1.H. Niknejad, H. Peirovi, M Jorjani, et al, "Properties of the Amniotic Membrane for Potential Use in Tissue Engineering", Eur Cell Mater, vol. 15, pp 88-99, 2008.
- 2.A. Roy, M. Mantay, C. Brannan, and S. Griffiths, "Placental Tissues as Biomaterials in Regenerative Medicine" BioMed Research International, vol. 2022, pp. 1-26, 2022.

**NVISION™**  
BIOMEDICAL TECHNOLOGIES

#### NVISION™

Biomedical Technologies  
MAR-ACT-001 Rev. B Nvision  
Biomedical Technologies, Inc.

#### ADDRESS

4590 Lockhill Selma Rd. San Antonio, TX 78249  
210.545.3713  
www.nvisionbiomed.com  
contact@nvisionbiomed.com



Innovative Wound Consultants

21222 Gathering Oak, Suite 103, San Antonio, TX 78260  
210.932.8038  
FAX: 210.579.6647  
www.innovativewoundconsultants.com



Trademarks and Registered Trademarks of Nvision Biomedical Technologies, Inc.  
Copyright ©2024 Nvision Biomedical Technologies, Inc. All rights reserved